

# OPERATING INSTRUCTIONS



## 510 Meta inverted confocal



**You must not operate this equipment without prior training from a BALM facility staff member.**

To arrange training and for help please contact:

Dr Ann Wheeler      ext: 2406      a.p.wheeler@qmul.ac.uk

Isma Ali              ext: 2407      i.ali@qmul.ac.uk

### Standard Operating Procedure

#### How to turn the equipment on:

1. Switch on the metal halide lamp
2. Switch on the Zeiss remote control
3. Switch on the computer and log in

#### How to turn the equipment off:

1. Switch off lasers and wait to cool
2. Switch off Zeiss remote control once lasers have cooled
3. Switch off computer
4. Switch off metal halide lamp

#### Rules of use:

This microscope should be treated with respect and care at all times.

This Microscope can only be used by Masters by Research or PhD students, Postdocs and members of staff.

The microscope lenses must be cleaned after every usage and the equipment treated carefully at all times.

If you have any problems at all with the microscope, no matter how trivial they may seem please see a technician immediately.

REMEMBER: You have 5GB of disk space on this microscope. Check before you start if you have room for your experiment. If not, delete your old data.

1) Open **ZEN** software



2) Click **Start System**




3) Place slide on microscope stage

4) In the **Ocular** tab  
choose an objective  
and a filter

5) Click **Online**  
and open the shutter

6) Find your sample



Use this button to change the objectives

Objectives available:


10x, 20x – air  
40x, 63x, 100x – oil

Automated, do not move by hand

**NB: please LOWER THE STAGE before changing between objectives (to avoid crashing lens onto slide)**

Use this button to change between the fluorescent filters

- 7) In the **Acquisition** tab select a **Configuration** and accept the prompt asking if you wish to switch on the lasers




- 8) In the **Laser** menu Switch the Argon laser from **standby** to **on**
- 9) In the **laser properties** (for the Argon laser) Set the **tube current** to 6.1A by adjusting the **output (%)**




NB. Ignore the warning about exceeding 50% output

- 10) Click **Auto Exposure** to generate an initial image

- 11) Click **Live** to see the live image, and select the **SPLIT** view to see the individual channels



- 12) Optimise the system by clicking **Merged**



Blue = Underexposed pixels

Red = Overexposed pixels

Adjust the **Gain** and **Digital Offset** so that the signal is within the dynamic range of the detectors (see below)

13) In the **Channels** menu click **Show all** and **Select all**

14) Set the pinhole to **1AU**


15) Alter the **Digital offset** so background is black not blue

16) Alter the **Gain (Master)** so there are no red pixels


17) Repeat for each channel



18) **Stop** the Live image aquisition





Click the **Merged** button again to return to the original colours



You can **Crop** your image and change the position, size or rotation of the cropped area by moving the crop box

Start scanning again and your cropped image will appear



19) In the **Acquisition Mode** menu set the image capture parameters


Recommended acquisition:

Frame Size: **1024 x 1024**


Speed: **7**

Averaging number: **4**


Bit depth: **8 bit**




20) Click **Snap**




21) Save your image on the **E: drive**  
(displayed in the top right of the workspace)




To export your images in TIF format, go to **File**, and **Export**.



Export as **Raw data – single plane** for overlay image.



Export as **Full resolution image window** for contents of the image window (will vary depending on which tab is open).



When you have finished, transfer all your data to the **Z network drive**

### PLEASE TIDY UP!!

Clean lenses, throw away used tissue/lens tissue, dispose of old slides in the yellow sharps bin